







Six good reasons to use aScope™ 3

-  aScope™ 3 offers improved patient safety as it is always available when needed
-  aScope™ 3 is sterile with no risk of cross-contamination
-  aScope™ 3 is portable, intuitive and fast to set-up and use
-  aScope™ 3 gives you 3 sizes in one system at no additional cost
-  aScope™ 3 is cost effective. Low investment costs. No repair and reprocessing costs
-  aScope™ 3 gives you easy access to practice your bronchoscopy skills



Ambu® aView™ Portable, intuitive monitor for Ambu® aScope™ 3 Family

Ambu® aView™ is a portable high resolution monitor that is compatible with all of Ambu's single-use flexible scopes in the Ambu® aScope™ 3 Family.

The clear and high resolution image quality enables easy navigation and fast identification of anatomical landmarks.

Key benefits of aView™ include


- Portable touch-screen monitor
- Operational in seconds
- Intuitive user interface
- Easy file transfer by USB
- Minimum 3 hours battery time
- Enhanced viewing with Image Extend
- Display of live image on external monitor, e.g. for training



Image Extend enlarges the live image on aView™ for better viewing at a distance.



Video-out functionality allows you to display the live image of aView™ on an external monitor, e.g. for training.

		SVU Placement of DLTs and BBs	Difficult intubation	Bronchial inspection	BAU/BW	Protected suction brush sampling	PDT	Management of retained secretion	Management of tough retained secretion	
	Ambu® aScope™ 3 Slim 3.8/1.2 Ref # 402001000	●	●	●						ETT ≥ 5 DLTs ≥ 37
	Ambu® aScope™ 3 Regular 5.0/2.2 Ref # 403001000	●	●	●	●	●	●	●		ETT ≥ 6 DLTs ≥ 39
	Ambu® aScope™ 3 Large 5.8/2.8 Ref # 414001000		●	●	●	●	●	●	●	ETT ≥ 7
	Ambu® aView™ Ref # 405002000									

The aScope™ 3 Family consists of 3 bronchoscope sizes. The colour coding makes it easy and fast to find the right size.

See for yourself at visualisation.ambu.com

Ambu A/S
Baltorpbakken 13
DK-2750 Ballerup
Denmark

Tel.: 1 410 768 6464
Fax: 1 800 262 8673
www.ambuUSA.com



496 3 605 01 V01 2016/03



Ambu® aScope™ 3 Family Tried and trusted

Ambu® aScope™ 3 Family ensures flexible bronchoscopes are always readily available



Better patient safety

aScope™ 3 ensures immediate access to a sterile flexible scope and eliminates the risk of cross-contamination, so you and your patients can feel safe.

Improved productivity

The single-use aScope™ 3 helps you save time and work smarter by eliminating time consuming steps required today to use, maintain and handle a reusable scope.

One system with three bronchoscope sizes to cover your needs

Whether you need a bronchoscope for intubation, placement of a DLT or suctioning of heavy secretion, the aScope™ 3 Family covers your needs. The colour coding makes it easy and fast to find the right size.



With the aScope™ 3 you always have a sterile bronchoscope available whenever you need it

aScope™ 3 is the safe and cost effective choice for fibreoptic intubation in the OR and ER



aScope™ 3 is ideal for a wide range of procedures in the OR:

- Intubation
- Airway inspection
- Placing and checking of Double Lumen Tubes
- Placement of Bronchial Blockers

Improving patient safety



You and your patients can feel safe as the aScope™ 3 ensures you always have a flexible scope readily available whenever needed.

Procedures made simple


The sterile, single-use and ready-to-go concept of aScope™ 3 simplifies flexible intubation. The aScope™ 3 system is highly portable, intuitive and fast to set up. After use there is no handling or reprocessing; simply discard the scope.

Cost effective


aScope™ 3 offers low investment costs and no repair and reprocessing costs. And you can stop worrying about surgery delays due to lack of an available flexible scope.



A clear image makes it easy to confirm that the DLT is correctly placed



Confirmation of patency through bronchial lumen of DLT



RUL apical segment
Maneuverability allows for easy access to the most difficult parts of the lung

aScope™ 3 makes lack of available bronchoscopes a thing of the past in the ICU



aScope™ 3 is ideal for a wide range of bedside procedures in the ICU:

- Airway inspection
- Difficult intubation
- Management of retained secretion
- Bronchial Aveolar Lavage
- Percutaneous Dilatational Tracheostomy
- Protected Specimen Brush sampling

Designed for bronchoscopy in ICU

Whether you need to intubate, perform a PDT or BAL procedure, or require high suction capacity for tough retained secretion, the three scope sizes in the aScope™ 3 Family cover your needs. The unique colour coding of the scopes makes it easy to select the right size. aScope™ 3 is immediately available, portable, and easy to use, making it the ideal choice for bedside bronchoscopy.

Your patient is safe

With the sterile single-use concept of aScope™ 3 your patient is safe as there is no risk of cross-contamination. And the ready-to-go concept ensures that you always have a scope readily at hand when you need it.

Cost effective

You can treat your patient here and now without having to wait for a scope to be found and brought to you. And you no longer have to worry about your scopes being broken and expensive repair costs.



Clear image quality makes aScope™ 3 ideal for bronchial inspection



aScope™ 3 can be navigated into all segments of the lung



aScope™ 3 used for Protected Specimen Brush sampling